


How to guide for the Juice Lab


Hydration Juices are supplied in cases of 12 x 0.5 litre cartons. Each 0.5 litres makes 10 litres of final product once mixed with water. Ambient storage of up to 12 months on unopened cartons. Open cartons have a 30 day refrigerated shelf life. When 2.5 litre fill line jug filled it has a 5 day refrigerated shelf life. When the Yorkshire dispenser is filled with juice it has a 3 day shelf life if refrigerated out of service times.



1 Unpack the graphics (back and base panels). With the 'man' to the left as you look at it gently fold out the side panels.



2 Place the base panel on the counter you are using and with the back panel upright align it on top. Ensure the upright panel sits on the base panel and not in the groove and the base folds up on the outside.



3 Holding the upright graphic peel the tape away from one section of the base. Once unpeeled the tape is very strong so make sure it is properly in position and lined up when you push it against the upright graphic.



4 Repeat to stick all three sides to the upright, position and add the copper riser. Make sure the riser is the front of the base unit with the drip tray completely flat on the counter.



5 Unpack dispensers, clean and rinse inside and out. Remove any packaging which may have settled in the tap. Ensure tap is in closed position. Separate and remove diffuser tube from ice tube. Fill ice tube with ice.



6 Place dispensers on riser. Take the right flavour labels, carefully remove one layer of the velcro dot, remove the tape from the reverse and position above the dispenser. Add the flavour label.



7 Use the 2.5L Fill Line Jug to mix the hydration juice. Fill the concentrate to 1+19 marker point, then fill the jug with water to the 2.5L marker point. Pour into dispenser, add the ice core and screw the metal lid on.



8 Repeat the filling of each dispenser and place on the riser lining up with the flavour labels. Add cups, information on pricing or how to use the Juice Lab as part of a meal deal.